



EARTH

POWERED BY LIGHT • INSPIRED BY SCIENCE

Solar light sparked life on Earth, and is the primary source of energy on our planet. Using the power of light for healing the body dates back to the Greeks, who observed that sunlight had a healing effect on humans. The Greeks called this form of therapy **heliotherapy**. They knew the benefits of using light to heal but they did not know that sunlight contains UV rays that if used consistently, are harmful to people and damage the tissue.

The Pioneer of Light Therapy

Dr. Niels Ryberg Finsen (1860–1904) delved into the theory of light as therapy, demonstrating the beneficial effects of various light wavelengths in the treatment of tuberculosis. Finsen was awarded the Nobel Prize in Medicine and Physiology in 1903. The use of concentrated light radiation for the treatment of smallpox, measles and lupus vulgaris opened a new avenue for medical science.


This same technology of light therapy has proven to be effective in the fields of dermatology, neurology and physiotherapy.



From the Nobel Prize to NASA

LED light therapy and the use of specific light waves for therapy has continued to advance. NASA noticed the increase in plant growth under specific light wavelengths. Since then, NASA began experimenting with LEDs to speed up the healing of wounds for pain management. The test results were positive and the research branched into the field of skin cancer treatment and cosmetics.

LED light therapy has been validated by thousands of studies because of its diversity and efficacy in the treatment of medical ailments.



Who We Are

Earth concept store will take you to the future where the natural power of light reverses the signs of aging and skin damage, making your home a dermatologist's clinic. The procedure is non-invasive, gentle, painless and safe; it takes no more than a few minutes a day at the comfort of your own home, at your convenience.

Our engineers have spent the last 5 years fine-tuning this technology—originally discovered by NASA—for our purposes, before even sending it to production, to be absolutely sure we can offer you the best product and results. All of our products are science-based, FDA approved and dermatologist recommended.

In order to prove the power and effectiveness of all Zero Gravity products, independent clinical studies were conducted, most notably and recently by Princeton Research Division in Princeton, New Jersey. The Princeton clinical trial was conducted with live subjects for a period of two months.

Our 5 Elements

Combining our 5 most important elements to provide comprehensive treatment with our latest technology.

Photo-
Biostimulation

Pain Relief
Management

Anti-Bacterial

Hair Growth

Special
Analysis

Photo-Biostimulation

While plants use light for photosynthesis, humans have a similar process called Photo Bio-stimulation, in which cells produce ATP (adenosine triphosphate), the fuel that powers life at a cellular level.

LED red light is known to promote production of cells, collagen fibers and elastin to rejuvenate, lift and tighten your skin. The red light is known to reduce melanin levels in the skin and improve the appearance of age spots and dark spots.

Perfectio Plus - Perfectio Plus will rejuvenate your skin appearance and structure using the dual action techniques – red LED light and topical heating care. A combination of powerful infrared LEDs offers deeper and more effective treatment to all skin layers.

Perfectio - Perfectio is a user-friendly product for acquiring younger, firmer skin at no time at all. Zero Gravity has perfected this revolutionary LED anti-aging facial rejuvenation technique, bringing it to the palm of your hand.



U.S. Patent No.
D 791,960
D 29/650,726

Perfectio X, Zero Gravity's flagship; a state of the art red and infrared LED light therapy device that provides incredible anti-aging benefits as well as treating chronic pain and decreasing inflammation. Its revolutionary glass probe is designed to use with cream or serum and to provide the best results.

Photo-Biostimulation – Perfectio X delivers precise dosage of energy to the skin to simultaneously produce collagen fibers and elastin, while gently erasing signs of aging in all skin layers. Naturally fills and erases fine lines and wrinkles, lifts and tightens the skin.

Pain Relief Management – The Perfectio X's powerful LED red and infrared beams penetrate deeply into the skin, muscles and joints. The light photons are absorbed by photoacceptors in your cells' mitochondria, which convert the light into ATP. This boost of energy increases the circulation in the blood vessels, micro-repairs tissues, and decreases inflammation and pain.



FDA
Cleared as a
Medical Device

CE
No.: A1860126 01
No.: A1860126 02

DERMATOLOGIST
RECOMMENDED
DERMATOLOGIST

APPROVED
PATENTED
APPROVED

CLINICALLY PROVEN
CLINICALLY PROVEN

D 29/657,936
D 16/013,068
D 29/665,452
ZL201830284928.6
ZL201830285765.3
2018-F-00633572

Anti - Bacterial

Sapphire's blue light therapy is a safe and effective treatment that uses dual action technology—blue light and topical heating care—to disinfect, detoxify and eliminate bacteria under the skin, thus engendering smoother, clearer and healthier skin.

Blue light therapy has also been proven to be extremely effective in the treatment of acute or chronic acne in all ages.

In clinical studies, more than 86% of study participants experienced a 74% reduction in acne.



FDA
Cleared as a
Medical Device

CE
No.TB170214760



U.S. Patent No.
D 605,762

Sapphire X blue light anti-bacterial device was developed with one main goal in mind: delivering the safest and most effective anti-bacterial blue light skin therapy treatment available.

The special waterproof glass head delivers an optimal light exposure for more effective treatment.

It offers a revolutionary approach to maintaining a healthy skin, using blue light therapy and topical heat to disinfect, detoxify and eliminate bacteria under the skin, for a smoother healthier complexion.



FDA
Cleared as a
Medical Device

CE
No.TB170214760



D 29/657,936
D 16/013,068
D 29/665,452
ZL201830284928.6
ZL201830285765.3
2018-F-00633572

Hair Growth

Research shows that 90% of men and 45% of women will suffer from loss of hair volume and density during their lives. In men, the process can lead to complete baldness, and in women to different levels of hair thinning in different parts of the scalp.

Grandio is a self-treatment device designed to treat hair loss, thinning hair and baldness. Grandio is suitable for both men and women, and is an effective method to boost hair regrowth. It provides a combination of three treatment mechanisms that effectively stop the balding process, neutralize the cause of rapid hair loss and encourage new hair growth.

Recreo is a hair growth device that uses low-level laser light therapy to reverse the gradual decline of hair production mechanisms—be the cause genetic or otherwise—and initiates a process of renewal and rehabilitation of the skin and hair follicles. It is unisex and it was recently approved by FDA as a medical device class II.



FDA
Cleared as a
Medical Device



U.S. Patent No.
D 29/657,967



CE

Earth's Luxury Accessories

Aecor Neck Devices are innovative beauty appliances that improve the overall appearance of the neck area. Various modes of magnetic heat lift and firm up the neck, promoting synthesis of protein and collagen and accelerating blood circulation for smoother look and feel. They utilize a non-invasive approach to micro-current technology that offers visibly apparent long lasting results.

Aecor Eye Device improves the overall appearance of the eyes and lips area. It utilizes a non-invasive approach to micro-current technology that offers visibly apparent long lasting results.




CE RoHS

Zero Gravity Collection

FDA
Cleared as a
Medical Device

CE
No.TB170214760



EARTH - ANTI-AGING AND THE SCIENCE BEHIND THE POWER OF LIGHT

Store Design Suggestion

EARTH - ANTI-AGING AND THE SCIENCE BEHIND THE POWER OF LIGHT







EARTH

EARTH - ANTI-AGING AND THE SCIENCE BEHIND THE POWER OF LIGHT



EARTH



EARTH